



# X2 GENERAL

EMC According to IEC 61326-1

SAFETY According to IEC 61010-1

X-RAY METERS STANDARD Complies with IEC 61674

**EXPOSURES NEEDED** One

**USB CABLES** 2 m (6.6 ft), 5 m (16.4 ft) and

5 m active extender

SIZE BASE UNIT

34 x 85 x 154 mm (1.3 x 3.3 x 6.1 in)

WEIGHT BASE UNIT

521 g (18.4 oz)

OPERATING TEMPERATURE
STORAGE TEMPERATURE

15 – 35 °C (59 – 95 °F) -25 – 70 °C (-13 – 158 °F)

POWER SOURCE

Rechargeable Li ion battery ~ 10 hours intensive usage

BATTERY TIME
BATTERY TESTED

According to UN 38.3

DISPLAY

4.3" LCD with capacitive touch ~ 10000 latest exposures

MEMORY SOFTWARE

X2 View for data handling and analysis. Also exports data to

Microsoft Excel.

PTB CERTIFICATE DE-17-M-PTB-0053

# X2 mAs

mAs

 RANGE
 0.001 – 9999 mAs

 RESOLUTION
 0.001 mAs

UNCERTAINTY

1%

mΑ

RANGE (PEAK)
RESOLUTION

0.1 – 1500 mA

UNCERTAINTY

0.01 mA 1%

TIME

**RANGE** 1 ms - 999 s **RESOLUTION** 0.1 ms

BANDWIDTH 1 kHz
UNCERTAINTY 0.5 %

**PULSES** 

**RANGE** 1 – 9999 pulses

**RESOLUTION** 1 pulse

**PULSE RATE** 

**RANGE** 0.1 – 200 pulses/s

**RESOLUTION** 0.1 pulse/s

mAs/PULSE

**RANGE** 0.001 – 9999 mAs

**RESOLUTION** 0.001 mAs

UNCERTAINTY 1%

WAVEFORM

**RESOLUTION** 125  $\mu$ s\* **BANDWIDTH** 1 kHz

#### UNFORS RAYSAFE UNCERTAINTY DEFINITION

The expanded uncertainty is stated as the combined uncertainty of measurement multiplied by the coverage factor k=2, which assuming a normal distribution has a coverage probability of 95 % (complies with GUM by ISO (1995, ISBN 92-67-10188-9)).

Instrument specifications are subject to purchased configuration. All specifications may change without notice.

<sup>\*</sup> automatically reduced for exposures longer than 3 s

## X2 R/F SFNSOR

**WEIGHT** 42 g (1.5 oz)

**SIZE** 14 x 22 x 79 mm (0.5 x 0.9 x 3.1 in)

**ACTIVE COMPENSATION** 

Beam quality independent for the following ranges:

 DOSE/DOSE RATE
 40 - 150 kVp, 1 - 14 mm Al HVL

 kVp
 40 - 150 kVp, up to 1 mm Cu

 TF
 60 - 120 kVp, up to 1 mm Cu

DOSE

 $\textbf{RANGE} \hspace{1.5cm} 1\,\text{nGy} - 9999\,\text{Gy}$ 

 $(0.1 \, \mu R - 9999 \, R)$ 

**UNCERTAINTY** 5 % or 5 nGy (0.5  $\mu$ R)

DOSE RATE

 $\textbf{RANGE} \hspace{1cm} 1 \, \text{nGy/s} - 500 \, \text{mGy/s}$ 

(5 μR/min – 3400 R/min)

 RESOLUTION
 1 nGy/s (5 μR/min)

 TRIG LEVEL
 50 nGy/s (340 μR/min)

UNCERTAINTY 5 % or

10 nGy/s (70 μR/min) x duty cycle

kVp

**RANGE** 40-150 kVp **MINIMUM DOSE**  $50 \text{ } \mu\text{Gy } (6 \text{ mR})$ 

**MINIMUM DOSE RATE (PEAK)** 10 μGy/s (70 mR/min)

UNCERTAINTY 2 %

HVL

**RANGE** 1-14 mm Al **MINIMUM DOSE**  $1 \mu \text{Gy} (120 \mu \text{R})$ 

MINIMUM DOSE RATE (PEAK)  $0.5 \mu Gy/s (3.5 mR/min) at > 70 kV$ 

2.5 μGy/s (17 mR/min) at 50 kV

UNCERTAINTY 10 %

TOTAL FILTRATION

RANGE 1.5 – 35 mm Al MINIMUM DOSE 50  $\mu$ Gy (6 mR) MINIMUM DOSE RATE (PEAK) 10  $\mu$ Gy/s (70 mR/min) UNCERTAINTY 10 % or 0.3 mm Al

TIME

 $\textbf{RANGE} \hspace{1cm} 1\,\text{ms} - 999\,\text{s}$ 

**RESOLUTION** 0.1 ms

**BANDWIDTH** 4 Hz – 4 kHz\*

UNCERTAINTY 0.5 %

\* automatically adjusted depending on signal level

**PULSES** 

**RANGE** 1 – 9999 pulses

**MINIMUM DOSE RATE (PEAK)**  $0.5 \mu Gy/s (3.5 mR/min)$ 

**PULSE RATE** 

**RANGE** 0.1 – 200 pulses/s

MINIMUM DOSE RATE (PEAK) 0.5 μGy/s (3.5 mR/min)

DOSE/PULSE

**RANGE** 1 nGy/pulse – 999 Gy/pulse

(0.1 μR/pulse – 999 R/pulse)

**MINIMUM DOSE RATE (PEAK)**  $0.5 \mu Gy/s (3.5 mR/min)$ 

WAVEFORMS

**RESOLUTION** 62.5  $\mu$ s\*

**BANDWIDTH kV**  $0.1 - 0.4 \text{ kHz}^{**}$ 

**BANDWIDTH DOSE RATE** 4 Hz – 4 kHz\*\*

\* automatically reduced for exposures longer than 1.5 s

<sup>\*\*</sup> automatically adjusted depending on signal level

## X2 MAM SENSOR

**WEIGHT** 42 g (1.5 oz)

**SIZE** 14 x 22 x 79 mm (0.5 x 0.9 x 3.1 in)

#### **ACTIVE COMPENSATION**

Beam quality independent for the following ranges:

#### DOSE/DOSE RATE & HVL

No selections needed.

With or without paddle, with or without phantom.

 Mo/Mo, Mo/Rh
 20 - 40 kVp

 Rh/Ag
 27 - 40 kVp

 Mo/Al, W/Rh, W/Ag,
 20 - 50 kVp

W/AI, Rh/Rh, Rh/AI

Mo/Cu, Rh/Cu, W/Cu, W/Ti 40 - 50 kVp

#### kVp

User selectable beam qualities.

Paddle compensation available when relevant.

**W/Ag** 20 - 40 kVp **W/Al** 20 - 50 kVp

Measuring above 40 kVp requires an X2 R/F Sensor + 2 mm Al (incl.)

**W/Rh** 20 – 40 kVp **Mo/Mo** 20 – 40 kVp

Mo/Rh 32 - 40 kVp using + 2 mm Al (incl)

**Rh/Ag** 27 - 40 kVp

Mo/Cu, W/Cu, W/Ti 40 - 50 kVp, using the X2 R/F Sensor

DOSE

**RANGE**  $1 \mu Gy - 9999 Gy (0.1 mR - 9999 R)$ 

UNCERTAINTY 5 %

DOSE RATE

**RANGE** 10  $\mu$ Gy/s - 300 mGy/s (70 mR/min - 2000 R/min)

UNCERTAINTY 5 %

kVp

**RANGE**  $20 - 50 \text{ kVp}^*$ 

Measuring above 40 kVp requires an X2 R/F Sensor

and on W/AI +2 mm AI (incl.)

MINIMUM DOSE 50 μGy (6 mR)

**MINIMUM DOSE RATE (PEAK)** 10 μGy/s (70 mR/min)

**UNCERTAINTY** 2 % or 0.5 kV (without paddle)

2 % or 0.7 kV (with paddle)

\* depending on beam quality, see active compensation

HVL

 $\begin{array}{lll} \mbox{RANGE} & 0.2 - 3.6 \ \mbox{mm Al} \\ \mbox{MINIMUM DOSE} & 1 \ \mbox{$\mu$Gy (0.1 \ \mbox{mR})$} \\ \mbox{UNCERTAINTY} & 5 \ \mbox{$\%$ above 25 kV} \\ & 10 \ \mbox{$\%$ below 25 kV} \\ \end{array}$ 

TIME

 RANGE
 1 ms - 999 s

 RESOLUTION
 0.1 ms

 BANDWIDTH
 400 Hz

 UNCERTAINTY
 0.5 %

**PULSES** 

**RANGE** 1 – 9999 pulses

**PULSE RATE** 

**RANGE** 0.1 - 200 pulses/s

DOSE/PULSE

**RANGE**  $1 \mu Gy/pulse - 999 Gy/pulse$ 

(0.1 mR/pulse - 999 R/pulse)

WAVEFORMS

 $\begin{array}{ll} \textbf{RESOLUTION} & 62.5~\mu\text{s}^{**} \\ \textbf{BANDWIDTH} & 400~\text{Hz} \end{array}$ 

 $<sup>^{**}</sup>$  automatically reduced for exposures longer than 1.5 s

# X2 LIGHT SENSOR

**WEIGHT** 136 g (4.8 oz)

**SIZE** 48 x 60 x 68 mm (1.9 x 2.4 x 2.7 in)

**CLASSIFICATION** DIN 5032 part 7 class B

**STANDARDS** Complies with relevant parts of

AAPM TG18, IEC 62563-1 and

IEC 61223-2-5.

LUMINANCE

**RANGE**  $0.01 - 10\ 000\ cd/m^2$ 

(0.03 - 34 000 fL)

**RESOLUTION** 0.001 cd/m<sup>2</sup> (0.001 fL)

APERTURE ANGLE 5°

MEASUREMENT AREA ø 10 mm (0.4 in)

UNCERTAINTY 3%

**ILLUMINANT A** 

**DEVIATION FROM** < 3 % (see figure Photopic Response)

HUMAN EYE V(λ) (f,')

ILLUMINANCE

**RANGE** 0.1 – 100 000 lux

(0.01 - 9000 fc)

**RESOLUTION** 0.01 lux (0.001 fc)

UNCERTAINTY 3%


**ILLUMINANT A** 

**DEVIATION FROM** < 3 % (see figure Photopic Response)


HUMAN EYE V(λ) (f,')

**COSINE DEVIATION (f<sub>2</sub>)** < 3 % (see figure Cosine Response)

### Photopic Response



### Cosine Response



## X2 SURVEY SENSOR

**DIMENSIONS** 14 x 66 x 192 mm (0.5 x 2.6 x 7.6 in)

**WEIGHT** 140 g (4.9 oz)

**ACTIVE COMPENSATION** H\*(10) – when selecting Sv

Air kerma – when selecting Gy or R

TRIG MODES

MANUAL Manual start and stop of measurement

AUTO Trig level (N80): 10  $\mu$ Gy/h (1.2 mR/h) or

20 μSv/h

H\*(10)

**RANGE** 0 nSv - 9999 Sv

**RESOLUTION** 1 nSv

**UNCERTAINTY** 10 %, N-series 20 – 150 kV

H\*(10) RATE

**RANGE**  $0 \mu Sv/h - 150 mSv/h$ 

**UNCERTAINTY** 10 % or 0.3  $\mu$ Sv/h, N-series 20 – 150 kV

AIR KERMA

**RANGE** 0 nGy – 9999 Gy (0 μR – 9999 R)

**RESOLUTION** 1 nGy  $(0.1 \mu R)$ 

**UNCERTAINTY** 5 %, RQA 50 – 150 kV

10 %. N-series 40 – 150 kV

AIR KERMA RATE

**RANGE**  $0 \mu Gy/h - 100 mGy/h (0 mR/h - 10 R/h)$ 

**UNCERTAINTY** 5 % or 0.3  $\mu$ Gy/h (0.03 mR/h),

RQA 50 - 150 kV

10% or 0.3  $\mu$ Gy/h (0.03 mR/h),

N-series 40 - 150 kV

**MEAN ENERGY** 

**RANGE** 30 – 120 keV

UNCERTAINTY 10 %

**MINIMUM DOSE RATE** 10  $\mu$ Sv/h or 10  $\mu$ Gy/h (1 mR/h)

**DEFINING STANDARD:** ISO 4037-1

TIME

**RANGE** 0.1 – 9999 s


**RESOLUTION** 0.01 s **BANDWIDTH** 1 Hz

WAVEFORM

**RESOLUTION** 10 ms **BANDWIDTH** 1 Hz

**MINIMUM DOSE RATE**  $1 \mu Sv/h \text{ or } 1 \mu Gy/h \text{ (0.1 mR/h)}$ 





# X2 CT SENSOR

WAVEFORMS RESOLUTION

BANDWIDTH

WEIGHT 86 g (3.0 oz) SIZE 14 x 22 x 219 mm (0.5 x 0.9 x 8.6 in) SIZE Ø 12.0 mm (0.47 in) **STANDARD** For measurements in accordance with IEC 60601-2-44 **ACTIVE LENGTH** 100 mm (3.94 in) < 5 % for 70 - 150 kV **ENERGY DEPENDENCE** (RQR, RQA and RQT beam qualities)  $55 - 110 \text{ kPa}, 15 - 35 ^{\circ}\text{C} (59 - 95 ^{\circ}\text{F})$ **AUTOMATIC ENVIRON-**MENTAL COMPENSATION DOSE RANGE  $10 \mu Gy - 999 Gy (1 mR - 999 R)$ 5 % UNCERTAINTY DOSE LENGTH PRODUCT RANGE 100 μGycm – 9999 Gycm (10 mRcm - 9999 Rcm) UNCERTAINTY 5 % **DOSE RATE** RANGE  $10 \mu Gy/s - 250 mGy/s$ (70 mR/min – 1700 R/min) 5 % UNCERTAINTY TIME RANGE 10 ms - 999 s RESOLUTION 1 ms 10 Hz **BANDWIDTH** 0.5 % UNCERTAINTY

1 ms

10 Hz

Unfors RaySafe offers comprehensive solutions for the X-ray room to measure the performance of X-ray equipment and to monitor medical staff dose in real-time.

RaySafe helps you avoid unnecessary radiation.

